


---

# Datasheet

## SERVO-n

---

GETINGE GROUP



# Contents

## SERVO-n

---

|  |    |                                      |    |
|--|----|--------------------------------------|----|
| General                                    | 3  | Special functions                    | 10 |
| The ventilator—General                     | 3  | Disconnection/Suction                | 10 |
| Ventilation—General                        | 4  | Monitoring and trends                | 11 |
| User interface                             | 4  | Alarms                               | 12 |
| Power supply                               | 5  | Autoset (alarm limits) specification | 12 |
| Gas supply                                 | 5  | Optional equipment                   | 13 |
| Operating conditions                       | 6  | Y sensor (option)                    | 14 |
| Non operating conditions                   | 6  | CO <sub>2</sub> analyzer (option)    | 15 |
| Standards—Safety and Functionality         | 6  | Edi module (option)                  | 16 |
| Display                                    | 7  | Aerogen nebulizers                   | 16 |
| Ventilation modes—Invasive ventilation     | 8  | Communication/Interface              | 17 |
| Invasive ventilation                       | 8  | Log function                         | 18 |
| Ventilation modes—Non invasive ventilation | 8  | Saving of data                       | 18 |
| Non invasive ventilation                   | 8  | Compressor Mini (option)             | 18 |
| High Flow therapy (option)                 | 9  | Service                              | 18 |
| Parameter settings                         | 9  | Ordering information                 | 18 |
| Backup parameter settings                  | 10 | Dimensional drawings                 | 19 |

# Technical specifications

## SERVO-n

### General

|                      |   |
|----------------------|---|
| Intended use         | <p>The SERVO-n ventilator system is:</p> <ul style="list-style-type: none"> <li>intended for respiratory support, monitoring and treatment of neonatal and pediatric patients</li> <li>to be used only by healthcare providers</li> <li>to be used only in professional healthcare facilities and for transport within these facilities</li> </ul> <p>For NAVA and Edi monitoring, it is in addition intended:</p> <ul style="list-style-type: none"> <li>to provide monitoring of the patient's breathing drive</li> <li>to improve synchrony between the ventilator system and patient when the electrical signal from the brain to the diaphragm is active</li> <li>for use on all patients with no contraindication for insertion/exchange of a nasogastric tube</li> </ul> |
| Instructions for use | Please carefully read the user's manual   |
| Legal manufacturer   | Maquet Critical Care AB   |
| Other products       | See separate datasheets.<br>Contact your local Getinge supplier for more information.   |

| The ventilator—General   | SERVO-n   | SERVO-n on mobile cart  |
|--|---|---|
| Base system weight   | Approximately 23 kg (50.7 lbs) <ul style="list-style-type: none"> <li>Patient unit 15 kg (33.0 lbs)</li> <li>User interface 4 kg (8.8 lbs)</li> <li>Handle 3 kg (6.6 lbs)</li> <li>Cable holder and cable 1 kg (2.2 lbs)</li> </ul> | Approximately 35 kg (77.2 lbs) <ul style="list-style-type: none"> <li>Base system approx. 23 kg (50.7 lbs)</li> <li>Mobile cart 12.5 kg (27.6 lbs)</li> </ul> |
| Dimensions of base (W x D), see dimensional drawings on page 19. | 368 x 205 mm (14.5" x 8.1")   | 647 x 547 mm (25.5" x 21.5")<br>incl. wheels  |
| Height (incl. user interface)                                    | 826 mm (32.5")  | 1368 mm (53.8")   |
| Wheels   | N/A   | Four wheels with separate brakes  |
| A-weighted sound pressure level ( $L_{PA}$ )                     | <40 dB, measured at a distance of 1 m (3.3 ft)  |   |
| A-weighted sound power level ( $L_{WA}$ )                        | <51 dB  |   |

# Technical specifications

## SERVO-n

### Ventilation—General

|  |   |
|--|---|
| Patient range                          | <ul style="list-style-type: none"> <li>• Neonatal: 0.3–8 kg (0.7–17.6 lbs)</li> <li>• Pediatric: 3–30 kg (6.6–66.0 lbs)</li> </ul>  |
| Bias flow                              | 0.5 l/min   |
| Internal compressible factor           | Max. 0.1 ml/cmH <sub>2</sub> O  |
| Gas delivery system                    | Microprocessor controlled valves  |
| Maximum airway pressure                | 125 cmH <sub>2</sub> O  |
| Method of triggering                   | Flow, pressure and Edi (with Edi module and Edi catheter)   |
| Inspiratory flow range                 | 0 to 33 l/min   |
| Pressure drop                          | <ul style="list-style-type: none"> <li>• Max. 6 cmH<sub>2</sub>O at a flow of 60 l/min (insp. channel)</li> <li>• Max. 3 cmH<sub>2</sub>O at a flow of 60 l/min (exp. channel)</li> </ul> |
| PEEP regulation                        | Microprocessor controlled valve   |
| Rise time, expiratory flow measurement | <12 ms for 10–90 % response at flow of 3–192 l/min  |
| Expiratory flow range                  | 0 to 192 l/min  |

### User interface

|              |  |
|--------------|--|
| Type         | TFT-LCD touchscreen                                    |
| Size         | 366 x 300 x 50 mm (14.4" x 11.8" x 2.0")               |
| Viewing area | 15" XGA, 1024 x 768 pixels with a 24-bit color palette |
| Weight       | Approximately 4 kg (8.8 lbs)                           |

# Technical specifications

## SERVO-n

### Power supply

Power supply, automatic range selection

- 100–120 V AC, 2 A, 50–60 Hz
- 220–240 V AC, 1 A, 50–60 Hz

Plug-in battery module

Battery backup (nickel-metal hydride, NiMH)

Six battery module slots. Two batteries are delivered with the ventilator.

Battery capacity

Rechargeable, 12 V, 3.5 Ah each

Battery backup time

Ranging from 60 minutes (2 batteries) to 180 minutes (6 batteries)

Recharge time

Approximately 3 h/battery

External 12 V DC

12.0 V–15.0 V DC, 10 A

Typical min. power consumption (no optional modules, no ongoing battery charging, normal panel backlight)

100 VA, 40 W at 230 V or 75 VA, 40 W at 110 V

Typical max. power consumption (with CO<sub>2</sub>), Edi and Y sensor modules, ongoing battery charging, max. panel backlight)

200 VA, 80 W at 230 V or 170 VA, 80 W at 110 V

### Gas supply

Inlet gas pressure air/O<sub>2</sub>

200–600 kPa / 2.0–6.0 bar / 29–87 PSI

Connection standards available

AGA, DISS, NIST, or French standard

Unavailable gas/loss of gas pressure

The flow from an unavailable gas (air or O<sub>2</sub>) is automatically compensated for so that the patient gets the preset volume and pressure.

Patient system gas connectors

Male 22 mm / female 15 mm. In accordance with ISO 5356-1

Gas exhaust port

Male 30 mm cone

# Technical specifications

## SERVO-n


### Operating conditions

|                                    |                            |
|------------------------------------|----------------------------|
| Temperature                        | 10 to 40 °C (50 to 104 °F) |
| Relative humidity                  | 15 to 95% non-condensing   |
| Atmospheric pressure               | 660 to 1060 hPa            |
| Lowest pressure in patient circuit | -400 cmH <sub>2</sub> O    |

### Non operating conditions

|                                    |                               |
|------------------------------------|-------------------------------|
| Temperature                        | -25 to +60 °C (-13 to 140 °F) |
| Relative humidity                  | < 95% condensing              |
| Atmospheric pressure               | 470 to 1060 hPa               |
| Lowest pressure in patient circuit | N/A                           |

### Standards— Safety and Functionality

|  |   |
|--|---|
| <br><small>0123</small> | The device complies with requirements and classification IIb of Medical Device Directive 93/42/EEC. CE Mark Notified Body number: 0123.   |
| Classification   | IEC 60601-1: 2005, Class I, continuous operation  |
| Standards  | <ul style="list-style-type: none"> <li>• ISO 80601-2-12:2011, ISO 80601-2-55:2011, EN 13544-1:2007+A1:2009</li> <li>• IEC 60601-1, Type B (equipment making physical contact with the patient and the gas pathways).</li> <li>• IEC 60601-1, Type BF (CO<sub>2</sub> analyzer, Y sensor, nebulizer patient unit and cable)</li> <li>• IEC 60601-1, Type CF-defibrillation proof (Edi catheter and cable)</li> </ul> |
| Ingress protection   | IP 21   |
| Electromagnetic compatibility (EMC)  | <ul style="list-style-type: none"> <li>• According to limits specified in IEC 60601-1-2: 2007.</li> <li>• The '<i>EMC Declaration, Information to the Responsible Organization</i>' is available from Maquet</li> </ul>   |

# Technical specifications

## SERVO-n

### Display

#### Views

- Basic view
- Advanced view
- Loops view
- Distance view
- Family view
- SERVO COMPASS® (option)

Each of the screen layout views offers a specific combination of displayed waveforms, loops and presented values.

#### Real-time waveforms

- Pressure
- Flow
- Volume
- CO<sub>2</sub> (with CO<sub>2</sub> analyzer option)
- Edi (with Edi module and Edi catheter)

#### Loops

- Pressure—Volume
- Pressure—Flow
- Volume—Flow

A reference loop and three overlaying loops can be displayed.

#### SERVO COMPASS

Visualizes volume (VT/BW) and pressure (total or driving) in relation to set targets in invasive modes.

# Technical specifications

## SERVO-n

### Ventilation modes—Invasive ventilation

|  |  |
|--|--|
| Controlled ventilation   | <ul style="list-style-type: none"> <li>• PC (Pressure Control)</li> <li>• VC (Volume Control)</li> <li>• PRVC (Pressure Regulated Volume Control)</li> </ul>   |
| Supported ventilation  | <ul style="list-style-type: none"> <li>• PS/CPAP (Pressure Support / Continuous Positive Airway Pressure)</li> <li>• VS (Volume Support)</li> </ul>  |
| Automode (option)  | <ul style="list-style-type: none"> <li>• Control mode: VC &lt;—&gt; Support mode: VS</li> <li>• Control mode: PC &lt;—&gt; Support mode: PS</li> <li>• Control mode: PRVC &lt;—&gt; Support mode: VS</li> </ul>                          |
| Combined ventilation   | <ul style="list-style-type: none"> <li>• SIMV (VC) + PS (Synchronized Intermittent Mandatory Ventilation)</li> <li>• SIMV (PC) + PS</li> <li>• SIMV (PRVC) + PS</li> <li>• Bi-Vent/APRV (Airway Pressure Release Ventilation)</li> </ul> |
| NAVA   | <ul style="list-style-type: none"> <li>• Neurally Adjusted Ventilatory Assist via endotracheal tube or tracheostomy</li> </ul>   |
| VC and SIMV (VC) + PS and Automode VC <—> VS are not available in the neonatal patient category. |  |

### Invasive ventilation

|                                 |          |          |
|---------------------------------|----------|----------|
| Max. leakage compensation level | Neonatal | 25 l/min |
|---------------------------------|----------|----------|

### Ventilation modes—Non invasive ventilation

|  |  |
|--|--|
| Controlled ventilation                                   | <ul style="list-style-type: none"> <li>• NIV PC (option)</li> </ul>  |
| Supported ventilation                                    | <ul style="list-style-type: none"> <li>• NIV PS (option)</li> <li>• Nasal CPAP</li> </ul>                                    |
| NIV NAVA   | <ul style="list-style-type: none"> <li>• Neurally Adjusted Ventilatory Assist via non-invasive patient interfaces</li> </ul> |
| NIV PS is not available in the neonatal patient category |  |

### Non invasive ventilation

|                                   |                                   |   |
|-----------------------------------|-----------------------------------|---|
| Max. leakage compensation level   | Pediatric and neonatal            | <ul style="list-style-type: none"> <li>• Inspiratory: up to 33 l/min</li> <li>• Expiratory: up to 25 l/min</li> <li>• Nasal CPAP: up to 20 l/min</li> </ul> |
| Disconnection flow (configurable) | Low<br>High<br>Disabled           | 7.5 l/min<br>15 l/min<br>Deactivates disconnection detection  |
| Connection detection              | Manual or automatic via bias flow |   |



# Technical specifications

## SERVO-n

### High Flow therapy (option)

|                    |                       |                              |
|--------------------|-----------------------|------------------------------|
| Flow setting range | Pediatric<br>Neonatal | 0.5–30 l/min<br>0.5–20 l/min |
|--------------------|-----------------------|------------------------------|

| Parameter settings                                    | Neonatal range | Pediatric range |
|---|----------------|-----------------|
| Tidal volume (ml)                                     | 2–50           | 10–350          |
| Minute volume (l/min)                                 | 0.1–7.5        | 0.3–20          |
| Apnea, time to alarm (s)                              | 1–45           | 2–45            |
| Max. apnea time in Automode (s)                       | 3–15           | 3–15            |
| Pressure level above PEEP (cmH <sub>2</sub> O)        | 0–80           | 0–80            |
| Pressure level above PEEP in NIV (cmH <sub>2</sub> O) | 0–60           | 0–60            |
| PEEP (cmH <sub>2</sub> O)                             | 1–50           | 1–50            |
| PEEP in NIV (cmH <sub>2</sub> O)                      | 2–20           | 2–20            |
| CPAP pressure (cmH <sub>2</sub> O)                    | 2–20           | 2–20            |
| Respiratory rate (breaths/min)                        | 4–150          | 4–150           |
| SIMV rate (breaths/min)                               | 1–60           | 1–60            |
| Breath cycle time, SIMV (s)                           | 0.5–15         | 0.5–15          |
| P <sub>High</sub> (cmH <sub>2</sub> O)                | 2–50           | 2–50            |
| T <sub>High</sub> (s)                                 | 0.2–30         | 0.2–30          |
| T <sub>PEEP</sub> (s)                                 | 0.1–10         | 0.1–10          |
| PS above P <sub>High</sub> (cmH <sub>2</sub> O)       | 0–78           | 0–78            |
| O <sub>2</sub> concentration (%)                      | 21–100         | 21–100          |
| I:E ratio   | 1:10–4:1       | 1:10–4:1        |
| Ti (s)  | 0.1–5          | 0.1–5           |
| NAVA level (cmH <sub>2</sub> O/μV)                    | 0–15           | 0–15            |
| Edi trigger (μV)                                      | 0.1–2.0        | 0.1–2.0         |
| T <sub>Pause</sub> (s)                                | –              | 0–1.5           |
| T <sub>Pause</sub> (% of breath cycle time)           | –              | 0–30            |
| Flow trigger (l/min)                                  | 0–0.5          | 0–0.5           |
| Pressure trigger (cmH <sub>2</sub> O)                 | –1 to –20      | –1 to –20       |
| Insp. rise time (% of breath cycle time)              | 0–20           | 0–20            |
| Insp. rise time (s)                                   | 0–0.2          | 0–0.2           |
| End inspiration (% of peak flow)                      | 1–70           | 1–70            |
| End inspiration (% of peak flow) in NIV               | 10–70          | 10–70           |

# Technical specifications

## SERVO-n

| Backup parameter settings                                    | Neonatal range | Pediatric range |
|--|----------------|-----------------|
| Inspiratory tidal volume (ml)                                | 2–50           | 10–350          |
| Pressure level above PEEP in backup (cmH <sub>2</sub> O)     | 5–80           | 5–80            |
| Pressure level above PEEP in NIV backup (cmH <sub>2</sub> O) | 5–60           | 5–60            |
| Respiratory rate in backup (breaths/min)                     | 4–150          | 4–150           |
| I:E ratio  | 1:10–4:1       | 1:10–4:1        |
| Ti (s)   | 0.1–5          | 0.1–5           |

| Special functions             | Setting range  |
|-------------------------------|--|
| Manual breath                 | Initiation of 1 breath (In SIMV mode initiation of 1 mandatory breath)     |
| Static measurements           | Insp. or exp. hold (0–30 seconds)  |
| Nebulization                  | 5–30 min/Continuous/Off  |
| O <sub>2</sub> boost level    | Off, 1–79 %  |
| O <sub>2</sub> boost function | Activate O <sub>2</sub> boost up to 1 minute                               |
| Leakage compensation          | On/Off   |
| Circuit compensation          | On/Off   |
| Edi monitoring                | In all ventilation modes and in Standby (with Edi module and Edi catheter) |
| Previous mode                 | Activates previously used mode   |
| Backup ventilation            | Backup On/Off  |
| Apnea management              | Several parameters   |

| Disconnection/Suction   |   |
|-------------------------|---|
| Pre-oxygenation time    | Max. 2 min                                |
| Post-oxygenation time   | Max. 1 min                                |
| Patient disconnected    | High priority alarm activated after 1 min |
| Adjustable oxygen level | 21–100 %                                  |

# Technical specifications

## SERVO-n

| Monitoring and trends                        |                                      |
|--|--------------------------------------|
| Peak airway pressure                         | Ppeak                                |
| Pause airway pressure                        | Pplat                                |
| Mean airway pressure                         | Pmean                                |
| Driving pressure                             | Pdrive                               |
| Positive end expiratory pressure             | PEEP                                 |
| Continuous positive airway pressure          | CPAP                                 |
| Spontaneous breaths per minute               | RR sp                                |
| Respiratory rate                             | RR                                   |
| Spontaneous expiratory minute volume         | MVe sp                               |
| Inspired minute volume                       | MVi                                  |
| Expired minute volume                        | MVe                                  |
| Leakage fraction (%)                         | Leakage                              |
| Inspired tidal volume                        | VTi                                  |
| Expired tidal volume                         | VTe                                  |
| End expiratory flow                          | Flowee                               |
| Measured oxygen concentration                | O <sub>2</sub> conc                  |
| CO <sub>2</sub> end tidal concentration      | etCO <sub>2</sub>                    |
| CO <sub>2</sub> minute elimination           | VCO <sub>2</sub>                     |
| CO <sub>2</sub> tidal elimination            | VTCO <sub>2</sub>                    |
| Dynamic compliance                           | Cdyn                                 |
| Static compliance                            | Cstatic                              |
| Inspiratory resistance                       | Ri                                   |
| Expiratory resistance                        | Re                                   |
| Work of breathing, ventilator                | WOBvent                              |
| Work of breathing, patient                   | WOBpat                               |
| Elastance                                    | E                                    |
| P 0.1  | P 0.1                                |
| Shallow Breathing Index                      | SBI                                  |
| Peak Edi value                               | Edipeak                              |
| Average Edipeak                              | Edipeak average (monitoring only)    |
| Average Edimin                               | Edimin average (monitoring only)     |
| Minimum Edi value                            | Edimin                               |
| Ratio of expired tidal volume to body weight | VT/BW                                |
| Switches to backup per minute                | Backup $\Sigma$ (trended value only) |
| Time in backup in percent per minute         | Backup % (trended value only)        |

# Technical specifications

## SERVO-n

| Alarms  | Neonatal range   | Pediatric range              |
|---|--|------------------------------|
| Airway pressure (upper alarm limit)               | 16–90 cmH <sub>2</sub> O   | 16–90 cmH <sub>2</sub> O     |
| Airway pressure NIV (upper alarm limit)           | 16–70 cmH <sub>2</sub> O   | 16–70 cmH <sub>2</sub> O     |
| Respiratory rate (upper and lower alarm limits)   | 1–160 breaths/min  | 1–160 breaths/min            |
| Expired minute volume (upper alarm limit)         | 0.02–30 l/min  | 0.02–30 l/min                |
| Expired minute volume (lower alarm limit)         | 0.01–20 l/min  | 0.01–20 l/min                |
| End expiratory pressure (upper alarm limit)       | 1–55 cmH <sub>2</sub> O  | 1–55 cmH <sub>2</sub> O      |
| End expiratory pressure (lower alarm limit)       | Off, 1–47 cmH <sub>2</sub> O   | Off, 1–47 cmH <sub>2</sub> O |
| No patient effort (Apnea) alarm                   | 1–45 s<br>Automatic return to support mode on patient triggering   | 2–45 s                       |
| No consistent patient effort                      | Yes, described in User's manual  |                              |
| High continuous pressure                          | Yes, described in User's manual  |                              |
| O <sub>2</sub> concentration                      | Set value ±5 vol% or ≤18 vol%  |                              |
| Gas supply  | Below 200 kPa (2.0 bar/29 PSI), above 600 kPa (6.0 bar/87 PSI)   |                              |
| Battery   | <ul style="list-style-type: none"> <li>Limited battery capacity: 10 min.</li> <li>No battery capacity: less than 3 min.</li> <li>Low battery voltage.</li> </ul> |                              |
| End tidal CO <sub>2</sub> (upper and lower limit) | 0.5–20 %, 4–100 mmHg, 0.5–14 kPa   |                              |
| Leakage too high                                  | Yes, described in User's manual  |                              |
| Technical   | Yes, described in User's manual  |                              |

| Autoset (alarm limits) specification                        | Invasive ventilation, controlled modes only                                 |
|---|---|
| High airway pressure  | Mean peak pressure +10 cmH <sub>2</sub> O or at least 35 cmH <sub>2</sub> O |
| Inspiratory tidal volume too high                           | The greater of VTi +30 % or VTi +2 ml                                       |
| Expiratory minute volume (upper alarm limit)                | Mean expiratory minute volume +50 %   |
| Expiratory minute volume (lower alarm limit)                | Mean expiratory minute volume –50 %   |
| Respiratory rate (upper alarm limit)                        | Mean respiratory rate +40 %   |
| Respiratory rate (lower alarm limit)                        | Mean respiratory rate –40 %   |
| End expiratory pressure (upper alarm limit)                 | Mean end expiratory pressure +5 cmH <sub>2</sub> O                          |
| End expiratory pressure (lower alarm limit)                 | Mean end expiratory pressure –3 cmH <sub>2</sub> O                          |
| End tidal CO <sub>2</sub> concentration (upper alarm limit) | Mean end tidal CO <sub>2</sub> concentration +25 %                          |
| End tidal CO <sub>2</sub> concentration (lower alarm limit) | Mean end tidal CO <sub>2</sub> concentration –25 %                          |

# Technical specifications

## SERVO-n

| Optional equipment          | Weight                | Dimensions  | Maximum load                                |
|-----------------------------|-----------------------|---|---|
| Mobile cart                 | 12.5 kg<br>(27.6 lbs) | W 647 x L 547 x H 557 mm<br>(W 25.5" x L 21.5" x H 21.9")   | -   |
| Drawer for mobile cart      | 0.6 kg<br>(1.3 lbs)   | W 247 x L 118 x H 302 mm<br>(W 9.7" x L 4.6" x H 11.9")   | -   |
| Shelf base                  | 2.5 kg<br>(5.5 lbs)   | W 207 x L 302 x H 43 mm<br>(W 8.2" x L 4.6" x H 1.7")   | -   |
| Pendant/bed holder          | 3.2 kg<br>(7.1 lbs)   | W 302 x L 302 x H 393 mm<br>(W 11.9" x L 11.9" x H 15.5")   | -   |
| Humidifier holder           | 0.5 kg<br>(1.1 lbs)   | W 243 x L 38 x H 185 mm<br>(W 9.6" x L 1.5" x H 7.3")   | 5 kg (11.0 lbs)                             |
| Support arm 178             | 2.5 kg<br>(5.5 lbs)   | L 900 mm (35.4")  | 1–3 kg ( 2.2–6.6 lbs)<br>depending on angle |
| User interface holder       | 0.6 kg<br>(1.3 lbs)   | W 46 x L 90 x H 123 mm<br>(W 1.8" x L 3.5" x H 4.8")  | -   |
| Cable holder for handle     | 0.1 kg<br>(0.2 lbs)   | W 138 x L 92 x H 155 mm<br>(W 5.4" x L 3.6" x H 6.1")   | 5 kg (11.0 lbs)                             |
| Waterbag/IV pole            | 0.4 kg<br>(0.9 lbs)   | W 148 x L 26 x H 1007 mm<br>(W 5.8" x L 1.0" x H 39.6")   | 1.5 kg (3.3 lbs)                            |
| Gas cylinder restrainer kit | 0.5 kg<br>(1.1 lbs)   | Upper: W 104 x L 65 x H 48 mm<br>(W 4.1" x L 2.5" x H 1.9")<br>Lower: W 106 x L 162 x H 76 mm<br>(W 4.1" x L 6.4" x H 3.0") | Two 4.5-liter bottles                       |
| Y piece holder              |                       | W 26 x L 52 x H 46 mm<br>(W 1.0" x L 2.0" x H 1.8")   |   |

# Technical specifications

## SERVO-n

| Y sensor (option)       | Size  | Weight          |
|-------------------------|---|-----------------|
| Y sensor module         | W 154 x L 90 x H 21 mm<br>(W 6.1" x L 3.5" x H 0.8")  | 280 g (0.6 lbs) |
| Y sensor                | W 18 x L 50 x H 27 mm<br>(W 0.7" x L 2.0" x H 1.1")   | 11 g            |
| Connectors and cables   | <ul style="list-style-type: none"> <li>• 15 mm male and female conical connector on flow sensor according to ISO 5356-1</li> <li>• Pressure port on module, pressure line, 2.0 m (6.6 ft), phthalate free PVC</li> <li>• Flow sensor cable, 2.0 m (6.6 ft)</li> </ul> |                 |
| Sensor material         | <ul style="list-style-type: none"> <li>• Single use: PC, Polycarbonate</li> <li>• Reusable: PEI, Polyetherimide or PSF, Polysulfone</li> </ul>  |                 |
| Power source            | Powered by the ventilator system, 4.5 W during normal operation.  |                 |
| Measuring method        | Hot Wire Anemometer (HWA)   |                 |
| Parameters              | <ul style="list-style-type: none"> <li>• Airway pressure.</li> <li>• Airway flow.</li> <li>• Inspiratory and expiratory volumes.</li> <li>• Trigger and End inspiration</li> </ul>  |                 |
| Measuring range         | <ul style="list-style-type: none"> <li>• Flow: 0.12 to 32 l/min</li> <li>• Pressure: -40 to 120 cmH<sub>2</sub>O</li> </ul>   |                 |
| Y sensor resistance     | 10 cmH <sub>2</sub> O/l/s at 30 l/min   |                 |
| Dead space              | ≤1 ml   |                 |
| Pressure line connector | Gable mounted bulk head connector to fit tubing with an inner diameter of 3–4 mm (0.12–0.16").  |                 |

# Technical specifications

## SERVO-n

| CO <sub>2</sub> analyzer (option)    | Size   | Weight               |
|--------------------------------------|--|----------------------|
| CO <sub>2</sub> analyzer module      | W 154 x L 90 x H 43 mm<br>(W 6.1" x L 3.5" x H 1.7")   | 0.45 kg (1.0 lbs)    |
| Sensor (Capnostat 5)                 | 32.0 x 47.0 x 21.6 mm<br>(1.3" x 1.9" x 0.8")  | 20 g                 |
| Operating temperature                | 10 °C to 33 °C   |                      |
| Airway adapter                       |  | 10 g                 |
| Power source                         | Powered by the ventilator  |                      |
| Connectors and cables                | Sensor   | 2.8 m (9.2 ft) cable |
| Measuring method                     | Mainstream, dual-wavelength, non-dispersive infrared   |                      |
| Parameters                           | <ul style="list-style-type: none"> <li>• CO<sub>2</sub> end tidal concentration (etCO<sub>2</sub>)</li> <li>• CO<sub>2</sub> minute elimination (VCO<sub>2</sub>)</li> <li>• CO<sub>2</sub> tidal elimination (VTCO<sub>2</sub>)</li> </ul>    |                      |
| Measuring range                      | <ul style="list-style-type: none"> <li>• 0 to 100 mmHg CO<sub>2</sub> partial pressure</li> <li>• 0 to 13.3 kPa CO<sub>2</sub> partial pressure</li> <li>• 0 to 13.2 % CO<sub>2</sub> volume (at a barometric pressure of 1013 hPa)</li> </ul> |                      |
| System response time CO <sub>2</sub> | The total system response time of the CO <sub>2</sub> monitor when exposed first to air and then to a gas mix with 5.0 % CO <sub>2</sub> is <250 ms.   |                      |
| Warm-up time                         | 15 s to initial CO <sub>2</sub> indication<br>maximum 2 minutes to full specification  |                      |
| Oxygen concentration compensation    | Automatic. Values supplied from the ventilator system  |                      |
| Barometric pressure compensation     | Automatic. Values supplied from the ventilator system  |                      |
| Digitizing rate                      | 100 Hz   |                      |
| Airway adapter dead space            | Neonatal/pediatric   | <1 cm <sup>3</sup>   |

# Technical specifications

## SERVO-n

| Edi module (option) | Size   | Weight            |
|---------------------|--|-------------------|
| Edi module          | W 154 x L 90 x H 21 mm<br>(W 6.1" x L 3.5" x H 0.8")   | 0.25 kg (0.6 lbs) |
| Edi catheter cable  | 2.0 m (6.6 ft)   | -                 |
| Power source        | Powered by the ventilator  |                   |
| Power consumption   | <3 W during normal operation   |                   |
| Parameters          | <ul style="list-style-type: none"> <li>• Edi waveform</li> <li>• ECG leads waveforms</li> <li>• NAVA estimated pressure waveform (Pest)</li> </ul> |                   |

| Aerogen nebulizers | Pro   | Solo  |
|--------------------|---|---|
| Size               | W 50 x L 50 x H 45 mm<br>(W 2.0" x L 2.0" x H 1.8") | W 48 x L 25 x H 67 mm<br>(W 1.9" x L 1.0" x H 2.6") |
| Weight             | Approx. 25 g  | Approx. 14 g  |
| Particle size      | 1–5 µm mass median aerodynamic diameter (MMAD)      |   |
| Flow rate          | >0.2 (average: ~0.4) ml/min                         |   |
| Max. volume        | 10 ml   | 6 ml  |
| Residual volume    | <0.1 ml for 3 ml dose                               |   |
| Control cable      | 1.8 m (5.9 ft)                                      |   |



# Technical specifications

## SERVO-n

| Communication/Interface             |   |
|-------------------------------------|---|
| Serial ports                        | <ul style="list-style-type: none"> <li>• Isolated</li> <li>• Two RS-232C ports. For data communication via the Servo Communication Interface (SCI).</li> </ul>  |
| Servo Communication Interface (SCI) | A protocol for data communication with external devices   |
| Alarm output connection (option)    | <ul style="list-style-type: none"> <li>• Isolated</li> <li>• 4-pin modular connector for communication of all active alarms</li> <li>• Switching capability: Max. 40 V DC, max. 500 mA, max. 20 W.</li> </ul>   |
| Data transfer via USB port          | <ul style="list-style-type: none"> <li>• Non-isolated</li> <li>• For transfer of trends, logs, screenshots and recordings to a USB memory stick</li> </ul>  |
| VGA port                            | <ul style="list-style-type: none"> <li>• Non-isolated</li> <li>• External monitors should be isolated from the ventilator system. This can be achieved if the connected monitor is powered via a medical grade transformer.</li> <li>• The secondary monitor output is not to be relied on for alarms.</li> </ul> |
| Ethernet port                       | <ul style="list-style-type: none"> <li>• Isolated</li> <li>• The network connection (LAN) port is for service use, and should only be used by personnel trained and authorized by Maquet.</li> </ul>  |
| MSync, HL7 converter (optional)     | See separate datasheet.   |

# Technical specifications

## SERVO-n

### Log function

Event log

- Alarms
- Ventilator settings
- Apnea periods
- Immediate functions

Diagnostic log

- Technical information
- Test results
- Service records
- Software installation
- Configuration information

### Saving of data

Recording of current waveform and parameter values

30 seconds of data will be recorded (15 seconds before and 15 seconds after activation). Up to 40 recordings can be stored

Saving screenshots

Up to 40 screenshots can be stored

Export files

Recordings, screenshots and the event log can be saved together in an export file and transferred to a USB memory stick

### Compressor Mini (option)

See separate datasheet

### Service

Regular maintenance

Once every 12 months or at least after 5000 operating hours

### Ordering information


SERVO-n, ventilator system and accessories: See separate information in “System flowchart, Ventilation, SERVO-n”, art no 68 85 739.

# Technical specifications


## SERVO-n

### Dimensional drawings


SERVO-n on mobile cart



SERVO-n holder



SERVO-n (patient unit) on SERVO-n holder



SERVO-n may be pending regulatory approvals to be marketed in your country. Contact your Getinge representative for more information.

This document is intended to provide information to an international audience outside of the US.

---

## GETINGE GROUP

**Legal manufacturer:**

Maquet Critical Care AB  
Röntgenvägen 2  
SE-171 54 Solna, Sweden  
Phone: +46 (0)10 335 73 00

[www.getingegroup.com](http://www.getingegroup.com)

---

**MAQUET**  
GETINGE GROUP

**ARJOHUNTLEIGH**  
GETINGE GROUP

**GETINGE**  
GETINGE GROUP

Getinge Group is a leading global provider of innovative solutions for operating rooms, intensive-care units, hospital wards, sterilization departments, elderly care and for life science companies and institutions. With a genuine passion for life we build quality and safety into every system. Our unique value proposition mirrors the continuum of care, enhancing efficiency throughout the clinical pathway. Based on our first-hand experience and close partnerships, we are able to exceed expectations from customers – improving the every-day life for people, today and tomorrow.