

Product Description:

一、 Use:

Drying box is mainly used for baking, drying, heat treatment and heating. Industrial and agricultural production, scientific research, health care and other units can be used in the laboratory, but not suitable for volatile substances, flammable materials, so as not to cause an explosion.

二、 Working principle

This machine adopts the thin steel plate is made in the studio sample shelf, an insulation layer between the studio and the studio at the bottom of the box body casing and is provided with an air hole, cold air from the bottom into the box, discharged from the top row of holes, two sets of heating components distributed in the lower part of the studio. Use as needed,

The temperature by the temperature sensor and temperature control control. The constant temperature system, the fuse is installed on the left side of the box body and has a movable door device, and the maintenance is convenient.

Technical Parameter:

Model	DTK-DO-25	DTK-DO-35	DTK-DO-45
Internal dimension WxHxD	250*250*250	350*350*350	350*350*350
External dimensionWxHxD	460*490*450	570*650*520	450*450*500
Model	DTK-DO-55	DTK-DO-75	DTK-DO-100
Internal dimension WxHxD	550*550*450	600*750*500	100*1600*600
External dimensionWxHxD	730*840*670	880*930*720	1140*2150*800
Power supply	AC220V/50HZ		
Temperature range	RT+10-250℃		
Temperature control accuracy	±1℃		
power	2200W		
Weight	120KG		
Sample tray	3 Block (3 layer)		
Alarm range	1~9999min		